



METHOD 150.3 DETERMINATION OF pH IN DRINKING WATER

> ADEM Surface Water Meeting October 23, 2019

> > Lennette West Hach Company

AGENDA

Why test for pH
What is pH
How to test for pH
Documentation
Probe care







WHY TEST FOR pH

• <u>ADEM</u>

Division 7 Regs – Drinking Water Revised August 9, 2019



335-7-2-.02 Sampling and Analytical Requirements. Samples to be used to demonstrate compliance with these regulations must be collected using procedures, containers, and preservatives established by EPA. Analysis of such samples must be performed using approved EPA methodology and by a laboratory certified by EPA or this Department which has demonstrated the ability to analyze the specific contaminants at an acceptable detectable limit established by EPA. Turbidity, chlorine residual, and secondary standards may be analyzed by a certified operator using procedures established by EPA. Confirmation samples may be required after the detection of a contaminant or the submittal of results which is questionable.



WHY TEST FOR pH

• <u>ADEM</u>

-Chapter 335-7-10



Be Riaht

- Operation, Record Keeping, And Reports -Section 335-7-10-.03
- (a) At surface water or ground water under the influence of surface water treatment plants serving community and NTNC systems, the following tests shall be performed and recorded during plant operation:
 - 6. pH of the raw and finished water each shift
- (b) At groundwater treatment facilities serving community and NTNC systems, the following tests shall be performed and recorded at least daily:
 - **5**. pH of the finished water from each source or plant if a chemical is fed to adjust the pH or if the pH is below 7.0;

WHY TEST FOR pH

• <u>ADEM</u>





- Control of Lead and Copper
- 2. The Langelier Index of the water is between -1.0 to +2,


pH is one of the parameters used in the calculation for Langelier Index





WHAT IS pH

- It is roughly the negative of the <u>logarithm</u> to base 10 of the <u>molar concentration</u>, measured in units of <u>moles per liter, of hydrogen ions</u>. More precisely it is the negative of the logarithm to base 10 of the <u>activity</u> of the hydrogen ion. (Wikipedia)
 - How acidic or basic (alkaline) is a substance.



HOW TO TEST FOR pH

• EPA Method 150.3 :

-Determination of pH in Drinking Water

-EPA 815-B-17-001 February 2017



Method 150.3: Determination of pH in Drinking Water





EPA METHOD 150.3 - SUMMARY

Can use any type meter

 Bench-top (lab)
 Portable (field)
 Continuous (process)



- Report pH with a resolution of 0.1 pH unit
- Proper probe maintenance is required
- Sample must be stirred or swirled no bubbles



EPA METHOD 150.3 - SUMMARY

- Reference pH Buffers

 Generally use 4, 7, and 10
 Adhere to expiration dates (Lot #s)
 - -Change out solutions weekly
- Sample Collection



Collect approximately 100ml of sample
Analyze immediately after collection





EPA METHOD 150.3 – CALIBRATION/VERIFICATION BENCHTOP AND PORTABLE

- <u>Calibration</u>
 - -Minimum of 2 buffer points (i.e. 7.0 and 10.0)
 - –At least weekly <u>or</u> if verification is out
- Verification
 - –After Cal, analyze a <u>different</u> pH buffer (i.e. 4.0)
 - Daily
 - –Must be within +/- 0.1 pH units





EPA METHOD 150.3 – CALIBRATION/VERIFICATION CONTINUOUS USING A PROBE (REMOVABLE)

<u>Calibration</u>

-Minimum of 2 buffer points (i.e., 7.0 and 10.0) -At least weekly or if verification is out

Verification

-After Cal, analyze a <u>different</u> pH buffer (i.e., 4.0) -Must be within +/- 0.1 pH units





DOCUMENTATION OF pH

Section 9 Quality Control

– "The laboratory or water system is required to maintain performance records that define the quality of data generated."

– Report both pH (nearest 0.1 unit) and temperature (°C) at time of analysis





pH – WHAT IS GOOD

- Record mv from each Buffer
- Record Slope (-59 mv)



Buffer	r	וע	
4.0	160	to	180
7.0	20	to	-20
10.0	-160	to	-180



ph probe care

- Consumable item
- Clean occasionally
- Stored properly
 - Do not let go dry
 - Storage Solution







DOCUMENTATION

- If it's not written down it wasn't done!
- How to prove something was done and done right - years later.
 - Calibration tell the instrument
 - Verification ensure the instrument is correct



Don't call the Nobel Committee just yet: We forgot to calibrate the instruments before the experiment...





pH Calibration Standard Methods 4500 H⁺ B-2011 (WW) EPA 150.3 (DW) OR Name: Date / Time: Expiration mV reading pH Buffer Lot Number Date Acceptable Range 4.0 +180 to +160 7.0 +20 to -20 10.0 -160 to -180 Known pH Slope (mV and %) Temp ^oC Sample ID Time pH Result Notes Known Should be 6.8 +/- 0.3



DOCUMENTATION

- Parameter testing / Method
- Name of technician/operator
- Date & Time of Cal &/or Verification
- Lot #(s) & Expiration Date of Reagents used
- Acceptable value or range
- QC's
- Sample ID
- Sample Result







QUESTIONS / COMMENTS

Lennette West Hach Company Iwest@hach.com